

Ceruloplasmin Microplate Assay Kit

Catalog # AS0059

Detection and Quantification of Ceruloplasmin Activity in Serum, Plasma, Other biological fluids Samples.

This instruction must be read in its entirety before using this product.

For research use only, Not for use in diagnostic procedures.

Contact information:

Tel:+1 (301) 446-2499 Fax:+1 (301) 446-2413

Email:techcn@signalwayantibody.com Web:www.sabbiotech.com



I. INTRODUCTION	2
II. KIT COMPONENTS	3
III. MATERIALS REQUIRED BUT NOT PROVIDED	3
VI.SAMPLE PREPARATION	4
V. ASSAY PROCEDURE	5
VI. CALCULATION	6
VII. TECHNICAL SUPPORT	7
VIII NOTES	7



I. INTRODUCTION

Ceruloplasmin is a copper containing protein found primarily in the blood. It carries approximately 70% of the total copper present in the blood. Ceruloplasmin exhibits a weak oxidase activity, which is a more accurate method of assessing ceruloplasmin in serum than immunodiffusion or other non-enzymatic assays. Normal Ceruloplasmin levels are generally 1-4 μ M (15-60 mg/dl), equivalent to approximately 50-150 mU/ml. Elevated amounts of serum ceruloplasmin are found in pregnancy, cancer, rheumatoid arthritis and in several mental conditions such as Alzheimer's, schizophrenia and OCD. Abnormally low amounts are found in Wilson's and Menkes' Diseases and in several other rare conditions.

The assay is initiated with the enzymatic catalysis of the PPD by CP. The enzyme catalysed reaction products can be measured at a colorimetric readout at 540 nm.



II.KIT COMPONENTS

Component	Volume	Storage
96-Well Microplate	1 plate	
Reaction Buffer	10 ml x 1	4 °C
Substrate	Powderx 1	4 °C, keep in dark
Plate Adhesive Strips	3 Strips	
Technical Manual	1 Manual	

Note:

Enzyme: add 10 ml distilled waterto dissolve before use.

III. MATERIALS REQUIRED BUT NOT PROVIDED

- 1. Microplate reader to read absorbance at 540 nm
- 2. Distilled water
- 3. Pipettor
- 4. Pipette tips
- 5. Mortar
- 6. Centrifuge
- 7. Timer
- 8. Ice



IV. SAMPLE PREPARATION

1.For serum or plasma samples

Detect directly.



V. ASSAY PROCEDURE

Add following reagents into the microplate:

Reagent	Sample	Blank
Sample	10 μΙ	
Distilled water		10 μΙ
Reaction Buffer	90 μΙ	90 μΙ
Substrate	100 μΙ	100 μΙ

Mix, put it into the oven at 37 °Cfor 10 minutes, measured at 540 nm and record the absorbance.



VI. CALCULATION

Unit Definition: one unit is defined as the enzyme that will oxidize 1 μ mol of Substrate per minute.

1. According to the volume of serum or plasma

CP (U/mI)=(OD_{Sample} -OD_{Blank})/
$$(\varepsilon \times d) \times V_{Total} / V_{Sample} / T$$

= 0.352×(OD_{Sample} -OD_{Blank})

ε: molar extinction coefficient,9.46× 10³L/mol/cm;

d: the optical path of 96-Well microplate, 0.6 cm;

V_{Total}: the total volume of the enzymatic reaction, 0.2 ml;

V_{Sample}: the volume of sample, 0.01 ml;

T: the reaction time, 10 minutes.



VII. TECHNICAL SUPPORT

For troubleshooting, information or assistance, please go online to www.sabbiotech.cn or contact us at techcn@signalwayantibody.com

VIII. NOTES