

## ANK3 Antibody

Catalog No: #37413

Package Size: #37413-1 50ul #37413-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ANK3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ANK3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human ankyrin 3, node of Ranvier (ankyrin G)
Target Name	ANK3
Other Names	ANKYRIN-G
Accession No.	Swiss-Prot#: Q12955NCBI Gene ID: 288Gene Accssion: NP_001191332
Concentration	2.6mg/ml
Formulation	Rabbit IgG in pH7.3 PBS, 0.05% NaN <sub>3</sub> , 50% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using #37413 at dilution 1/50.

## Background

Ankyrins are a family of proteins that are believed to link the integral membrane proteins to the underlying spectrin-actin cytoskeleton and play key roles in activities such as cell motility, activation, proliferation, contact, and the maintenance of specialized membrane domains. Multiple isoforms of ankyrin with different affinities for various target proteins are expressed in a tissue-specific, developmentally regulated manner. Most ankyrins are

typically composed of three structural domains: an amino-terminal domain containing multiple ankyrin repeats; a central region with a highly conserved spectrin binding domain; and a carboxy-terminal regulatory domain which is the least conserved and subject to variation. Ankyrin 3 is an immunologically distinct gene product from ankyrins 1 and 2, and was originally found at the axonal initial segment and nodes of Ranvier of neurons in the central and peripheral nervous systems. Multiple transcript variants encoding different isoforms have been found for this gene.

## Published Papers

---

et al., Genetic Alterations in Papillary Thyroid Carcinoma With Hashimoto's Thyroiditis O ANK3, an Indolent Maintainer of Papillary Thyroid Carcinoma. In Front Oncol  
on 2022 May 12 by Chao Zeng, Jiali Long, et al.. PMID:35646694, , (2022)  
[PMID:35646694](#)

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.