

## TIPARP Antibody

Catalog No: #37797

Package Size: #37797-1 50ul #37797-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TIPARP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total TIPARP protein.
Immunogen Type	Peptide
Immunogen Description	Recombinant fusion protein of human TIPARP.
Target Name	TIPARP
Other Names	PARP7; ARTD14; pART14
Accession No.	Swiss-Prot#: Q7Z3E1NCBI Gene ID: 25976Gene Accssion: NP_056323
SDS-PAGE MW	76kd
Concentration	2.9mg/ml
Formulation	Rabbit IgG in pH7.3 PBS, 0.05% NaN <sub>3</sub> , 50% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:200-1:1000

## Images



Gel: 6%SDS-PAGE  
 Lysates (from left to right): Mouse spleen tissue  
 Amount of lysate: 40ug per lane  
 Primary antibody: 1/200 dilution  
 Secondary antibody dilution: 1/8000  
 Exposure time: 5 minutes

## Background

This gene encodes a member of the poly(ADP-ribose) polymerase superfamily. Studies of the mouse ortholog have shown that the encoded protein catalyzes histone poly(ADP-ribosyl)ation and may be involved in T-cell function. Alternative splicing results in multiple transcript variants.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.