

## Quiescin Q6 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C06382Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Quiescin Q6 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Quiescin Q6 QSOX1
Conjugates	Cy5
Target Name	Quiescin Q6
Other Names	hQSOX; Human quiescin Q6 mRNA, partial cds; Q6 antibody QSCN6; QSOX1; QSOX1_HUMAN; Quiescin Q6; Quiescin Q6 sulfhydryl oxidase 1; Skin sulfhydryl oxidase; Sox; Sulfhydryl oxidase 1; Sulfhydryl oxidase 1 precursor.
Accession No.	NCBI Gene ID:5768
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

IF:1:50-200

## Background

This gene encodes a protein that contains domains of thioredoxin and ERV1, members of two long-standing gene families. The gene expression is induced as fibroblasts begin to exit the proliferative cycle and enter quiescence, suggesting that this gene plays an important role in growth regulation. Two transcript variants encoding two different isoforms have been found for this gene.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.