

Histone H3(Phospho-Tyr41) Rabbit Polyclonal Antibody

Catalog No: #12088

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #12088-1 50ul #12088-2 100ul


Description

Product Name	Histone H3(Phospho-Tyr41) Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	The Histone H3(Phospho Tyr41) Rabbit Polyclonal Antibody detects endogenous Histone H3(Phospho Tyr41) protein.
Immunogen Type	peptide
Immunogen Description	Peptide sequence around phosphorylation site of tyrosine 41(H-R-Y(p)-R-P) derived from Human Histone H3.
Target Name	Histone H3(Phospho-Tyr41)
Modification	Phospho
Other Names	H3 histone antibody; HIST1H3A antibody; Histone cluster 1; H3a antibody
Accession No.	Swiss-Prot#:P68431
SDS-PAGE MW	15-17kd
Concentration	1.0mg/ml
Formulation	Rabbit Antiserum, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000~2000

Images



Western blot analysis of extracts from HeLa cells, untreated (-) or treated with Nocodazole (1µg/ml, 24 hr; +), using #12088 diluted at 1:2,000.

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.