

WDR16 Polyclonal Antibody Cy5.5 Conjugated

Catalog No: #C06583Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WDR16 Polyclonal Antibody Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 435-485 620 derived from human WDR16
Conjugates	Cy5.5
Target Name	WDR16
Other Names	WDR16; WDRPUH; WD repeat-containing protein 16; WD40-repeat protein up-regulated in HCC
Accession No.	Swiss-Prot#:Q8N1V2NCBI Gene ID:146845
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

IF:1:50-200

Background

WD repeat containing proteins, such as WDR16, play crucial roles in a wide range of physiologic functions, including signal transduction, RNA processing, remodeling the cytoskeleton, regulation of vesicular traffic, and cell division.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.