

# JAKMIP2 Polyclonal Antibody Cy5.5 Conjugated

Catalog No: #C06621Cy5.5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JAKMIP2 Polyclonal Antibody Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF(IHC-P)
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human JAKMIP2 NECC1
Conjugates	Cy5.5
Target Name	JAKMIP2
Other Names	CTCL tumor antigen HD CL 04; CTCL tumor antigen HD-CL-04; Jak and microtubule interacting protein 2; JAKMIP 2; JAKMIP2; JAMIP 2; Jamip2; Janus kinase and microtubule interacting protein 2; Janus kinase and microtubule-interacting protein 2; JKIP2; JKIP2_HUMAN; KIAA0555; NECC 1; NECC1; Neuroendocrine
Accession No.	NCBI Gene ID:84525
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

## Application Details

IF:1:50-200

## Background

The protein encoded by this gene is reported to be a component of the Golgi matrix. It may act as a golgin protein by negatively regulating transit of secretory cargo and by acting as a structural scaffold of the Golgi. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012].

Note: This product is for in vitro research use only and is not intended for use in humans or animals.