

COLEC10 Polyclonal Antibody Cy5.5 Conjugated

Catalog No: #C07176Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	COLEC10 Polyclonal Antibody Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 142-187 277 derived from human COLEC10
Conjugates	Cy5.5
Target Name	COLEC10
Other Names	CLL1; Collectin-10; Collectin liver protein 1; CL-L1; Collectin-34; CL-34; COLEC10; UNQ366 PRO702
Accession No.	Swiss-Prot#:Q9Y6Z7NCBI Gene ID:10584
Cell Localization	Secreted
Concentration	1mg ml
Formulation	Aqueous buffered solution containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Background

Lectin that binds to various sugars: galactose > mannose = fucose > N-acetylglucosamine > N-acetylgalactosamine.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.