

TNFRSF10D antibody

Catalog No: #38726

Package Size: #38726-1 50ul #38726-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TNFRSF10D antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total TNFRSF10D protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TNFRSF10D (NP_003831.2).
Target Name	TNFRSF10D
Other Names	TNFRSF10D;CD264;DCR2;TRAIL-R4;TRAILR4;TRUNDD
Accession No.	Uniprot:Q9UBN6GeneID:8793
SDS-PAGE MW	47KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of THP-1 cells, using TNFRSF10D antibody.



Immunohistochemistry of paraffin-embedded mouse kidney using TNFRSF10D Rabbit pAb.

Background

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular TRAIL-binding domain, a transmembrane domain, and a truncated cytoplasmic death domain. This receptor does not induce apoptosis, and has been shown to play an inhibitory role in TRAIL-induced cell apoptosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.