

DTYMK antibody

Catalog No: #38855

Package Size: #38855-1 50ul #38855-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	DTYMK antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total DTYMK protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human DTYMK.
Target Name	DTYMK
Other Names	CDC8; TMPK; TYMK; PP3731;
Accession No.	Swiss-Prot#: P23919NCBI Gene ID: 1841
SDS-PAGE MW	23kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details


Western blotting: □ 1:500 - 1:2000

Immunofluorescence: 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using DTYMK antibody at 1:1000 dilution.



Immunofluorescence analysis of A549 cells using DTYMK antibody. Blue: DAPI for nuclear staining.

Background

Catalyzes the conversion of dTMP to dTDP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.