

# NCR1 antibody

Catalog No: #38916

Package Size: #38916-1 50ul #38916-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NCR1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total NCR1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human NCR1.
Target Name	NCR1
Other Names	LY94; CD335; NKP46; NK-p46;
Accession No.	Swiss-Prot#: O76036NCBI Gene ID: 9437
SDS-PAGE MW	34kd
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □1:500 - 1:2000

Immunofluorescence: □1:50 - 1:100

## Background

Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.