

TPPP3 antibody

Catalog No: #39171

Package Size: #39171-1 50ul #39171-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TPPP3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total TPPP3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TPPP3 (NP_057224.2).
Target Name	TPPP3
Other Names	TPPP3;CGI-38;TPPP/p20;p20;p25gamma
Accession No.	Uniprot:Q9BW30GeneID:51673
SDS-PAGE MW	19KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details


WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse spinal cord, using TPPP3 antibody.

Western blot analysis of extracts of various cell lines, using TPPP3 antibody.



Background

Regulator of microtubule dynamic that has microtubule bundling activity. Required for embryo implantation; possibly by regulating beta-catenin (By similarity). Also required for decidualization via regulation of beta-catenin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.