

## VPS26 Antibody

Catalog No: #39420

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VPS26 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Human
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein
Target Name	VPS26
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

## Application Details

WB

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.