

# ZNF215 Antibody

Catalog No: #39610



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |   |
|-----------------------|---|
| Product Name          | ZNF215 Antibody   |
| Host Species          | Rabbit  |
| Purification          | Affinity purified   |
| Applications          | ELISA WB  |
| Species Reactivity    | Hu  |
| Immunogen Type        | Recombinant Protein   |
| Immunogen Description | Raised against a recombinant human ZNF215 protein 1-252aa (BC014538).   |
| Target Name           | ZNF215  |
| Concentration         | 1.0mg/ml  |
| Formulation           | PBS, pH 7.4 with 0.02% Sodium Azide   |
| Storage               | This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing.<br>Expiration date is one (1) year from date of receipt. |

## Application Details

ELISA WB

## Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.