

ZNF331 Antibody

Catalog No: #39637



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF331 Antibody
Host Species	Rabbit
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant Protein
Immunogen Description	Raised against a recombinant human ZNF331 protein 75-352aa (BC009433).
Target Name	ZNF331
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal testis showing nucleus staining with anti-ZNF331 antibody at a dilution of 1/100.



Western Blot analysis of fetal liver Lysate with anti-ZNF331 antibody

WB: 1:1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.