

ZFHX3 Antibody

Catalog No: #39695



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZFHX3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Rt Ms
Immunogen Type	Recombinant Protein
Immunogen Description	raised against a recombinant human ZFHX3 protein 517-787aa (BC029653).
Target Name	ZFHX3
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA WB IHC

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal brain showing nuclear staining with anti-ZFHX3 antibody at a dilution of 1/100.



Western Blot analysis of HeLa, MCF-7 and A549 cell Lysate with anti-ZFHX3 antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.