

ZNF33B Antibody

Catalog No: #39756



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ZNF33B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Specificity	Human and predicted against mouse and rat.
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal ZNF33B (1) antibody was raised against a recombinant human ZNF33B protein 302-556aa (BC109232).
Target Name	ZNF33B
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:1000-1:2000

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.