

LRRC8A Antibody

Catalog No: #39905

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRRC8A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB IHC
Species Reactivity	Hu Ms Rt
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal LRRC8A (1) antibody was raised against a recombinated human LRRC8A protein 648-803aa (BC051322).
Target Name	LRRC8A
Concentration	1.0mg/ml
Formulation	PBS, pH 7.4 with 0.02% Sodium Azide
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.


Application Details

ELISA: 1:20000-1:80000 WB: 1:500-1:1000 IHC: 1:100-1:500

Images



Immunohistochemical staining of formalin-fixed paraffin-embedded fetal kidney showing cytoplasmic staining with anti-LRRC8A antibody at a dilution of 1/100



Note: This product is for in vitro research use only and is not intended for use in humans or animals.