

Recombinant Human Angiotensin-Converting Enzyme 2 ACE2 (C-6His)

Catalog No: #AP75124

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human Angiotensin-Converting Enzyme 2 ACE2 (C-6His) |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q9BYF1 |
| Calculated MW | 84.63kDa |
| Target Sequence | <p>QSTIEEQAKTFLDKFNHEAEDLFYQSSLASWNYNTNITEENVQNMNAGDKWSAFLKEQSTLAQMYPLQEIQ NLTVKLQLQALQQNGSSVLSSEDKSKRLNTILNTMSTIYSTGKVCNPDNPQECLEPGLNEIMANSLDYNERLW AWESWRSEVGKQLRPLYEYVVLKNEMARANYEDYGDYWRGDYEVNGVDGYDYSRGLIEDVEHTFEEIK PLYEHLHAYVRAKLMNAYPSYISPIGCLPAHLLGDMWGRFWTNLYSLTVPFQKPNIDVTDAMVDQAWDAQR IFKEAEKFFVSVGLPNMTQGFWENSMLTDPGNVQKAVCHPTAWDLGKGDFRILMCTKVTMDDFLTAHHEMG HIQYDMAYAAQPFLLRNANEGFHEAVGEIMSLSAATPKHLKSIGLLSPDFQEDNETEINFLKQALTIVGTLPF TYMLEKWRWVMVFKGEIPKDQWMKKWEMKREIVGVPEVPHDETYCDPASLFHVSNDYSFIRYYTRTYQF QFQEALCQAAKHEGPLHKCDISNTEAGQKLFNMLRLGKSEPWTALENVGAKNMNVRPLLNYFEPLFTWL KDQNKNSFVGWSTDWSPYADQSIKVRISLKSALGDKAYEWNNDNEMYLFRSSVAYAMRQYFLKVKNQMILFGE EDVRVANLKRISFNFFVTAPKNVSDIIPRTEVEKAIRMSRSRINDAFRLNDNSLEFLGIQPTLGPPNQPPVSD HHHHHH</p> |
| Formulation | Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 0.1mM ZnCl ₂ , 10% Glycerol, pH 7.5. |
| Storage | Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.