

## VAV2 Antibody

Catalog No: #40401

Package Size: #40401-1 50ul #40401-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	VAV2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total VAV2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human vav 2 guanine nucleotide exchange factor
Target Name	VAV2
Other Names	VAV-2
Accession No.	Swiss-Prot:P52735 Gene Accssion:NP_003362
Concentration	2.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry:1:10-1:50

## Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #40401 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #40401 at dilution 1/20.

## Background

VAV2 is the second member of the VAV guanine nucleotide exchange factor family of oncogenes. Unlike VAV1, which is expressed exclusively in hematopoietic cells, VAV2 transcripts were found in most tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.?

Note: This product is for in vitro research use only and is not intended for use in humans or animals.