

VWCE Antibody

Catalog No: #40403



Package Size: #40403-1 50ul #40403-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VWCE Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total VWCE protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human von Willebrand factor C and EGF domains
Target Name	VWCE
Other Names	VWC1; URG11
Accession No.	Swiss-Prot:Q96DN2 Gene Accssion:NP_689931
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry:1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #40403 at dilution 1/30.



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #40403 at dilution 1/30.

Background

May be a regulatory element in the beta-catenin signaling pathway and a target for chemoprevention of hepatocellular carcinoma.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.