

Recombinant Human Netrin(C-6His)

Catalog No: #AP75402



Package Size: #AP75402-1 10ug #AP75402-2 50ug #AP75402-3 500ug #AP75402-4 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Netrin(C-6His)
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Calculated MW	43.5kDa
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.