

## HDAC5/9 (Phospho-Ser259/220) Polyclonal Antibody

Catalog No: #12226



Package Size: #12226-1 50ul #12226-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HDAC5/9 (Phospho-Ser259/220) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-HDAC5/9 (S259/220) Polyclonal Antibody detects endogenous levels of HDAC5/9 protein only when phosphorylated at S259/220.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC5/9 around the phosphorylation site of S259/220.
Target Name	HDAC5/9
Modification	Phospho
Other Names	HDAC5; KIAA0600; Histone deacetylase 5; HD5; Antigen NY-CO-9; HDAC9; HDAC7; HDAC7B; HDRP; KIAA0744; MITR; Histone deacetylase 9; HD9; Histone deacetylase 7B; HD7; HD7b; Histone deacetylase-related protein; MEF2-interacting transcription rep
Accession No.	Swiss-Prot: Q9UQL6/Q9UKV0NCBI Gene ID: 10014/9734
Target Species	human
SDS-PAGE MW	121kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300


ELISA: 1/10000

Not yet tested in other applications

## Images



Western Blot analysis of HepG2 cells using Phospho-HDAC5/9 (S259/220) Polyclonal Antibody



Western Blot analysis of COLO HELA 293T cells using Phospho-HDAC5/9 (S259/220) Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.