

Cdk9 (Phospho-Thr186) Polyclonal Antibody

Catalog No: #12242



Package Size: #12242-1 50ul #12242-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cdk9 (Phospho-Thr186) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Cdk9 (T186) Polyclonal Antibody detects endogenous levels of Cdk9 protein only when phosphorylated at T186.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Cdk9 around the phosphorylation site of T186.
Target Name	Cdk9
Modification	Phospho
Other Names	CDK9; CDC2L4; TAK; Cyclin-dependent kinase 9; C-2K; Cell division cycle 2-like protein kinase 4; Cell division protein kinase 9; Serine/threonine-protein kinase PITALRE; Tat-associated kinase complex catalytic subunit
Accession No.	Swiss-Prot: P50750NCBI Gene ID: 1025
Target Species	human
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1/500 - 1/2000

ELISA: 1/5000

Not yet tested in other applications

Images

Western Blot analysis of 3T3 cells using Phospho-Cdk9 (T186) Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.