

## hnRNP K (Phospho-Ser284) Polyclonal Antibody

Catalog No: #12243



Package Size: #12243-1 50ul #12243-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | hnRNP K (Phospho-Ser284) Polyclonal Antibody  |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.     |
| Applications          | WB IHC ELISA  |
| Species Reactivity    | Hu Ms Rt Mk   |
| Specificity           | Phospho-hnRNP K (S284) Polyclonal Antibody detects endogenous levels of hnRNP K protein only when phosphorylated at S284. |
| Immunogen Type        | peptide   |
| Immunogen Description | Synthesized peptide derived from human hnRNP K around the phosphorylation site of S284.                                   |
| Target Name           | hnRNP K   |
| Modification          | Phospho   |
| Other Names           | HNRNPK; HNRPK; Heterogeneous nuclear ribonucleoprotein K; hnRNP K; Transformation up-regulated nuclear protein; TUNP      |
| Accession No.         | Swiss-Prot: P61978NCBI Gene ID: 3190  |
| Target Species        | human   |
| SDS-PAGE MW           | 51kd  |
| Concentration         | 1mg/ml  |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| Storage               | Store at -20°C/1 year   |

## Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/20000

Not yet tested in other applications

## Images

Western Blot analysis of 3T3 cells using Phospho-hnRNP K (S284) Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.