

# Recombinant Human Stimulator of Interferon Genes Protein, Sting (N-Sumo-6His)

Catalog No: #AP76379

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Stimulator of Interferon Genes Protein, Sting (N-Sumo-6His)
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q86WV6
Calculated MW	33.8kDa
Target Sequence	MGHHHHHHGSLQDSEVNQEAKPEVKPEVKPETHINLKVSDGSSEIFFKIKKTTPLRRLMEAFKRQGGKEMDSL RFLYDGIRIQADQAPEDLDMEDNDIIEAHREQIGGVAHGLAWSYYIGYLRLLPELQARIRTYNQHYNNLLRGAV SQRLYILLPLDCGVPDNLSMADPNIRFLDKLPQQTGDHAGIKDRVYSNSIYELLENGQRAGTCVLEYATPLQTL FAMSQYSQAGFSREDRLEQAKLFCRTLEDILADAPESQNNCRLLIAYQEPADDSSFSLSQEVLRHLRQEEKEEV
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Store at < -70°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.