

## Syk (Phospho-Tyr525) Polyclonal Antibody

Catalog No: #12294



Package Size: #12294-1 50ul #12294-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Syk (Phospho-Tyr525) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	Phospho-Syk (Y525) Polyclonal Antibody detects endogenous levels of Syk protein only when phosphorylated at Y525.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Syk around the phosphorylation site of Y525.
Target Name	Syk
Modification	Phospho
Other Names	SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk
Accession No.	Swiss-Prot: P43405NCBI Gene ID: 6850
Target Species	human
SDS-PAGE MW	72kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1/500 - 1/2000


Immunohistochemistry: 1/100 - 1/300

ELISA: 1/20000

Not yet tested in other applications

## Images

Western Blot analysis of COS7 cells using Phospho-Syk (Y525) Polyclonal Antibody



---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.