

## CD71 (Phospho-Ser24) Polyclonal Antibody

Catalog No: #12325



Package Size: #12325-1 50ul #12325-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD71 (Phospho-Ser24) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms
Specificity	Phospho-CD71 (S24) Polyclonal Antibody detects endogenous levels of CD71 protein only when phosphorylated at S24.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human CD71 around the phosphorylation site of S24.
Target Name	CD71
Modification	Phospho
Other Names	TFRC; Transferrin receptor protein 1; TR; TfR; TfR1; Trfr; T9; p90; CD antigen CD71
Accession No.	Swiss-Prot: P02786NCBI Gene ID: 4155
Target Species	human
SDS-PAGE MW	89kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/10000

Not yet tested in other applications

## Images



Western Blot analysis of 293 cells using Phospho-CD71 (S24) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.