

## JAK3 (Phospho-Tyr785) Polyclonal Antibody

Catalog No: #12331



Package Size: #12331-1 50ul #12331-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	JAK3 (Phospho-Tyr785) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-JAK3 (Y785) Polyclonal Antibody detects endogenous levels of JAK3 protein only when phosphorylated at Y785.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human JAK3 around the phosphorylation site of Y785.
Target Name	JAK3
Modification	Phospho
Other Names	JAK3; Tyrosine-protein kinase JAK3; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK
Accession No.	Swiss-Prot: P52333NCBI Gene ID: 3718
Target Species	human
SDS-PAGE MW	125kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

ELISA: 1/5000

Not yet tested in other applications

## Images



Western Blot analysis of 453 cells using Phospho-JAK3 (Y785) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.