

Gab 1 (Phospho-Tyr659) Polyclonal Antibody

Catalog No: #12359



Package Size: #12359-1 50ul #12359-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Gab 1 (Phospho-Tyr659) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Gab 1 (Y659) Polyclonal Antibody detects endogenous levels of Gab 1 protein only when phosphorylated at Y659.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Gab 1 around the phosphorylation site of Y659.
Target Name	Gab 1
Modification	Phospho
Other Names	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1
Accession No.	Swiss-Prot: Q13480NCBI Gene ID: 2549
Target Species	human
SDS-PAGE MW	76kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300


Immunofluorescence: 1/200 - 1/1000

ELISA: 1/5000

Not yet tested in other applications

Images

Western Blot analysis of 3T3 cells using Phospho-Gab 1 (Y659) Polyclonal Antibody



Note: This product is for in vitro research use only and is not intended for use in humans or animals.