

Catenin- β (Phospho-Ser552) Polyclonal Antibody

Catalog No: #12395



Package Size: #12395-1 50ul #12395-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Catenin- β (Phospho-Ser552) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Phospho-Catenin- β (S552) Polyclonal Antibody detects endogenous levels of Catenin- β protein only when phosphorylated at S552.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Catenin- β around the phosphorylation site of S552.
Target Name	Catenin- β
Modification	Phospho
Other Names	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin
Accession No.	Swiss-Prot: P35222NCBI Gene ID: 1499
Target Species	human
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1/500 - 1/2000

ELISA: 1/40000

Not yet tested in other applications

Images



Western Blot analysis of NIH-3T3, cells using Phospho-Catenin- β (S552) Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.