

## Histone H2B (Acetyl-Lys126) Polyclonal Antibody

Catalog No: #HW116



Package Size: #HW116-1 50ul #HW116-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Histone H2B (Acetyl-Lys126) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-Histone H2B (K126) Polyclonal Antibody detects endogenous levels of Histone H2B protein only when acetylated at K126.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the human Histone H2B around the acetylation site of K126.
Target Name	Histone H2B
Modification	Acetyl
Other Names	HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF;
Accession No.	Swiss-Prot: Q96A08NCBI Gene ID: 255626
Target Species	human
SDS-PAGE MW	15kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western blotting: 1/500 - 1/2000


ELISA: 1/20000

Not yet tested in other applications

## Images



Western Blot analysis of A549 cells using Acetyl-Histone H2B (K126) Polyclonal Antibody



Western Blot analysis of A549 cells using Acetyl-Histone H2B (K126) Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.