

TRADD antibody

Catalog No: #22562

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRADD antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu Ms
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TRADD
Target Name	TRADD
Accession No.	NCBI Gene ID: 8717NCBI mRNA#: NM_003789NCBI Protein#: NP_003780
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 34kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:200

Images



Western blot analysis on COS7 and HuvEc cell lysate using TRADD Antibody

Background

The protein encoded by this gene is a death domain containing adaptor molecule that interacts with TNFRSF1A/TNFR1 and mediates programmed cell death signaling and NF-kappaB activation. This protein binds adaptor protein TRAF2, reduces the recruitment of inhibitor-of-apoptosis proteins (IAPs) by TRAF2, and thus suppresses TRAF2 mediated apoptosis. This protein can also interact with receptor TNFRSF6/FAS and adaptor protein FADD/MORT1, and is involved in the Fas-induced cell death pathway. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.