

## GAP-43 Monoclonal Antibody

Catalog No: #40485



Package Size: #40485-1 50ul #40485-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	GAP-43 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous GAP-43 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Target Name	GAP-43
Accession No.	Swiss-Prot: P17677NCBI Gene ID: 2596
SDS-PAGE MW	38,43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western blotting: 1:1000-2000

Immunohistochemistry: 1:200-500

## Images



Western blot analysis of 1) HeLa, 2) 293T, 3) Mouse Brain Tissue, 4) Rat Brain Tissue using GAP-43 Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using GAP-43 Monoclonal Antibody.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.