

Rabbit IgG (Light Chain Specific) Monoclonal Antibody

Catalog No: #40498



Package Size: #40498-1 50ul #40498-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Rabbit IgG (Light Chain Specific) Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Rt
Specificity	The antibody reacts with the light chains on rabbit IgG and with those common to other rabbit immunoglobulins. The antibody does not react with the heavy chain of rabbit IgG. The antibody minimal cross-reaction with bovine, goat, Armenian hamster, horse,
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-5000

Images



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.