

Angptl7 Polyclonal Antibody

Catalog No: #40588



Package Size: #40588-1 50ul #40588-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angptl7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse
Specificity	Angptl7 Polyclonal Antibody detects endogenous levels of Angptl7 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Angptl7.
Target Name	Angptl7
Other Names	ANGPTL7; CDT6; Angiopoietin-related protein 7; Angiopoietin-like factor; Angiopoietin-like protein 7; Cornea-derived transcript 6 protein
Accession No.	Swiss-Prot: O43827NCBI Gene ID: 10218
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


WB 1:500 - 1:2000.

IHC 1:100 - 1:300.

ELISA: 1:10000..

IF 1:50-200

Images



Western Blot analysis of SH-SY5Y cells using Angptl7 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.