

# Vasopressin-neurophysin 2-copeptin

Catalog No: #AP78485

Package Size: #AP78485-1 50ug #AP78485-2 100ug #AP78485-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vasopressin-neurophysin 2-copeptin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	32-124AA
Other Names	ARVP,VP,AVP-NPII
Accession No.	P01185Gene name:AVP
Calculated MW	10.23
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Vasopressin has a direct antidiuretic action on the kidney, it also causes vasoconstriction of the peripheral vessels.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.