

Vitronectin

Catalog No: #AP78771



Package Size: #AP78771-1 50ug #AP78771-2 100ug #AP78771-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vitronectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	S-protein, Serum-spreading factor, V75
Accession No.	P04004 Gene name: VTN
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Somatomedin-B is a growth hormone-dependent serum factor with protease-inhibiting activity.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.