

EGFR Polyclonal Antibody

Catalog No: #40867



Package Size: #40867-1 50ul #40867-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EGFR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	EGFR Polyclonal Antibody detects endogenous levels of EGFR protein.
Immunogen Description	Synthesized peptide derived from human EGFR around the non-phosphorylation site of T693.
Target Name	EGFR
Other Names	EGFR; ERBB; ERBB1; HER1; Epidermal growth factor receptor; Proto-oncogene c-ErbB-1; Receptor tyrosine-protein kinase erbB-1
Accession No.	Swiss-Prot: P00533NCBI Gene ID: 1956
SDS-PAGE MW	175kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using EGFR Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.