

IL-8 Polyclonal Antibody

Catalog No: #41063



Package Size: #41063-1 50ul #41063-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Hu
Specificity	IL-8 Polyclonal Antibody detects endogenous levels of IL-8 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-8.
Target Name	IL-8
Other Names	IL8; CXCL8; Interleukin-8; IL-8; C-X-C motif chemokine 8; Emotakin; Granulocyte chemotactic protein 1; GCP-1; Monocyte-derived neutrophil chemotactic factor; MDNCF; Monocyte-derived neutrophil-activating peptide; MONAP; Neutrophil-activati
Accession No.	Swiss-Prot: P10145NCBI Gene ID: 3576
SDS-PAGE MW	8kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


IF 1:50-200

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using IL-8 Polyclonal Antibody

Published Papers

el at., Mechanisms of regenerative preconditioning in pigs with subtotal hepatectomies. In Cir Cir on 2022 by Orue-Echebarria, Javier Vaquero, et al..PMID:35944117, , (2022)

[PMID:35944117](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.