

JAK3 Polyclonal Antibody

Catalog No: #41082



Package Size: #41082-1 50ul #41082-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	JAK3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu Ms Rt
Specificity	JAK3 Polyclonal Antibody detects endogenous levels of JAK3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human JAK3 around the non-phosphorylation site of Y785.
Target Name	JAK3
Other Names	JAK3; Tyrosine-protein kinase JAK3; Janus kinase 3; JAK-3; Leukocyte janus kinase; L-JAK
Accession No.	Swiss-Prot: P52333NCBI Gene ID: 3718
SDS-PAGE MW	125kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using JAK3 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.