

Lamin B1 Polyclonal Antibody

Catalog No: #41100



Package Size: #41100-1 50ul #41100-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lamin B1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	Lamin B1 Polyclonal Antibody detects endogenous levels of Lamin B1 protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Lamin B1.
Target Name	Lamin B1
Other Names	LMNB1; LMN2; LMNB; Lamin-B1
Accession No.	Swiss-Prot: P20700NCBI Gene ID: 4001
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of mouse brain cells using Lamin B1 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.