

# NFκB-p100 Polyclonal Antibody

Catalog No: #41224



Package Size: #41224-1 50ul #41224-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NFκB-p100 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt
Specificity	NFκB-p100 Polyclonal Antibody detects endogenous levels of NFκB-p100 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human NFκB-p100 around the non-phosphorylation site of S869.
Target Name	NFκB-p100
Other Names	NFKB2; LYT10; Nuclear factor NF-kappa-B p100 subunit; DNA-binding factor KBF2; H2TF1; Lymphocyte translocation chromosome 10 protein; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2; Oncogene Lyt-10; Lyt10
Accession No.	Swiss-Prot: Q00653NCBI Gene ID: 4791
SDS-PAGE MW	97kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.


ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of Jurkat cells using NFκB-p100 Polyclonal Antibody



Western Blot analysis of Jurkat cells using NFκB-p100 Polyclonal Antibody

## Published Papers

el at., Single-cell analysis reveals the immune heterogeneity and interactions in lungs undergoing hepatic ischemia-reperfusion. *Int Immunopharmacol*. 2023 Nov;103:107111. doi: 10.1016/j.intimp.2023.107111. PMID: 37844464. (2023)  
[PMID:37844464](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.