

THAP12 Antibody

Catalog No: #47277

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	THAP12 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total THAP12 protein.
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human THAP12
Target Name	THAP12
Other Names	DAP4; THAP0; PRKRIR; P52rIPK
Accession No.	Swiss-Prot#:O43422NCBI Gene ID:5612Gene Accssion:BC021992
Concentration	0.8
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 30-150

Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47277(THAP12 Antibody) at dilution 1/45. (Original magnification: 200)

Background

Upstream regulator of interferon-induced serine/threonine protein kinase R (PKR). May block the PKR-inhibitory function of DNAJC3, resulting in restoration of kinase activity and suppression of cell growth.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.