

TIPRL Antibody

Catalog No: #47284

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TIPRL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms, Rt
Specificity	The antibody detects endogenous levels of total TIPRL protein.
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human TIPRL
Target Name	TIPRL
Other Names	TIP; TIP41; TIPRL1
Accession No.	Swiss-Prot#:O75663NCBI Gene ID:261726Gene Accssion:BC009506
Concentration	1.1
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

IHC dilution:1: 50-300

Images



The image is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using 47284(TIPRL Antibody) at dilution 1/65.(Original magnification: 200)



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using 47284(TIPRL Antibody) at dilution 1/65.(Original magnification: 200)

Background

TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.