

JPT2 Antibody

Catalog No: #47286

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	JPT2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total JPT2 protein.
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human JPT2
Target Name	JPT2
Other Names	L11; HN1L; C16orf34
Accession No.	Swiss-Prot#:Q9H910NCBI Gene ID:90861Gene Accssion:BC014438
Concentration	1.4
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 30-150

Images



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 47286(JPT2 Antibody) at dilution 1/80.(Original magnification: 200)

Background

The function of this protein remains unknown.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.