

## ZNF580 Antibody

Catalog No: #47346

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF580 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total ZNF580 protein.
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human ZNF580
Target Name	ZNF580
Accession No.	Swiss-Prot#:Q9UK33NCBI Gene ID:51157Gene Accssion:BC017698
Concentration	1.3
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

IHC dilution:1: 50-300

## Images



The image is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using 47346(ZNF580 Antibody) at dilution 1/75.(Original magnification: 200)



The image is immunohistochemistry of paraffin-embedded Human breast cancer tissue using 47346(ZNF580 Antibody) at dilution 1/75.(Original magnification: 200)

## Background

---

Involved in the regulation of endothelial cell proliferation and migration. Mediates H<sub>2</sub>O<sub>2</sub>-induced leukocyte chemotaxis by elevating interleukin-8 production and may play a role in inflammation. May be involved in transcriptional regulation.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.