

TMEM9 Antibody

Catalog No: #47353

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total TMEM9 protein.
Immunogen Type	Peptide
Immunogen Description	Fusion protein of human TMEM9
Target Name	TMEM9
Other Names	DERM4; TMEM9A
Accession No.	Swiss-Prot#:Q9P0T7NCBI Gene ID:252839Gene Accssion:BC001106
Concentration	1.1
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 50-300

Images



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 47353(TMEM9 Antibody) at dilution 1/50.(Original magnification: 200)

Background

May be involved in intracellular transport.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.