

NUDT8 Antibody

Catalog No: #47511

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUDT8 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu, Ms
Specificity	The antibody detects endogenous levels of total NUDT8 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human NUDT8
Target Name	NUDT8
Accession No.	Swiss-Prot#:Q8WV74NCBI Gene ID:254552Gene Accssion:NP_001230679
Concentration	0.7
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

IHC dilution:1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 47511(NUDT8 Antibody) at dilution 1/25.(Original magnification: 200)

Background

Probably mediates the hydrolysis of some nucleoside diphosphate derivatives.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.