

## C14orf180 Antibody

Catalog No: #47536

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	C14orf180 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total C14orf180 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide of human C14orf180
Target Name	C14orf180
Other Names	NRAC; C14orf77
Accession No.	Swiss-Prot#:Q8N912NCBI Gene ID:400258Gene Accssion:NP_001008404
Concentration	0.4
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

IHC dilution:1: 25-100

## Images



The image is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 47536(C14orf180 Antibody) at dilution 1/20.(Original magnification: 200)

## Background

C14orf180 is a multi-pass membrane protein. The function of the C14orf180 protein is not known.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.